Buffered Peptone Water

Diluent and non-selective pre-enrichment liquid medium for microbiological examination of food, according to ISO 6887, 11290, 21528 and 6579.

DESCRIPTION

Buffered Peptone Water (BPW) is a liquid medium recommended by ISO 6579 for increasing the recovery of injured *Salmonella* spp. from food and associated samples prior to selective enrichment and isolation.

According to ISO 21528, BPW is used for detection or enumeration of Enterobacteriaceae within foodstuffs.

Used as diluent, BPW complies with ISO 6887 and 11290 for the enumeration of organisms.

TYPICAL FORMULA

	(g/l)
Enzymatic Digest of Casein	10.0
Sodium Chloride	5.0
Disodium Hydrogen Phosphate	3.5*
Potassium Dihydrogen Phosphate	1.5
Final pH 7.0 ± 0.2 at 25°C	

*Equivalent to 9.0 g of disodium hydrogen phosphate dodecahydrate.

METHOD PRINCIPLE

Enzymatic digest of casein provides amino acids, nitrogen, carbon and minerals. Sodium chloride maintains the osmotic balance of the medium. Phosphates are the buffering agents.

PREPARATION

Dehydrated medium Suspend 20.0 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 min.

TEST PROCEDURE

Suspend the sample in BPW to make dilutions as required.

For pre-enrichment, add sample to BPW at a ratio of 1:10 or 1:9 depending on the method being used. Incubate at 37 ± 1°C for 16-20 hours before transfer to selective enrichment media.

INTERPRETING RESULTS

Turbidity indicates microbial growth.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.

Prepared medium: clear, light amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and tubes at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Prepared medium: 2 years.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC tables.

Inoculum for use as diluent: 10³-10⁴ CFU.

Incubation conditions: 45 min - 1 h / 18-27°C.

QC Table 1.

Microorganism	WDCM	Specification
<i>Escherichia coli</i>	WDCM 00012	± 30% colonies of original count
<i>Staphylococcus aureus</i>	WDCM 00034	± 30% colonies of original count
<i>Listeria monocytogenes</i> 4b	WDCM 00021	± 30% colonies of original count

Inoculum for productivity: ≤ 100 CFU.

Incubation conditions: 18 ± 2 h / 37 ± 1°C.

QC Table 2.

Microorganism	WDCM	Specification
<i>Salmonella typhimurium</i>	WDCM 00031	Good growth/turbidity of the medium
<i>Salmonella enteritidis</i>	WDCM 00030	Good growth/turbidity of the medium
<i>Escherichia coli</i>	WDCM 00012	Good growth/turbidity of the medium

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use and must be used only by properly trained operators.

DISPOSAL OF WASTE


Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 11290-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Listeria monocytogenes* and *Listeria* spp. – Part 2: Enumeration method.
3. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 1: Detection of *Enterobacteriaceae*.
4. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique.
5. ISO 6579-1:2017. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* spp. – Part 1: Detection of *Salmonella* spp.
6. Rose (2001) Isolation and identification of *Salmonella* from meat, poultry and egg products. In Microbiology laboratory guidebook, 3rd ed., Food Safety and Inspection Service, U.S. Department of Agriculture, Washington, D.C.
7. ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
8. Sadowski (1977) J. Food Technol. 12:85.
9. Edel and Kampelmacher (1973) Bull. W.H.O. 48:167.

PRESENTATION	Category	Packaging	Ref.
Buffered Peptone Water	Media in tubes	20 x 9 ml	24199
Buffered Peptone Water	Media in tubes	100 x 9 ml	26199
Buffered Peptone Water	Media in tubes	20 x 10 ml	24099
Buffered Peptone Water (Double Concentration)	Media in tubes	20 x 9 ml	24463
Buffered Peptone Water	Media in bottles	6 x 90 ml	414030
Buffered Peptone Water	Media in bottles	25 x 90 ml	454030
Buffered Peptone Water	Media in bottles	6 x 200 ml	412090
Buffered Peptone Water	Media in bottles	6 x 225 ml	414020
Buffered Peptone Water	Media in bottles	25 x 225 ml	451402
Buffered Peptone Water - Bags	Media in bags	3 x 3 l	499030
Buffered Peptone Water - Bags	Media in bags	3 x 5 l	499035
Buffered Peptone Water	Dehydrated media	100 g	621014
Buffered Peptone Water	Dehydrated media	500 g	611014
Buffered Peptone Water	Dehydrated media	5 kg	6110145

TABLE OF SYMBOLS

LOT Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 11290-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Listeria monocytogenes* and *Listeria* spp. – Part 2: Enumeration method.
3. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 1: Detection of *Enterobacteriaceae*.
4. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique.
5. ISO 6579-1:2017. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* spp. – Part 1: Detection of *Salmonella* spp.
6. Rose (2001) Isolation and identification of *Salmonella* from meat, poultry and egg products. In Microbiology laboratory guidebook, 3rd ed., Food Safety and Inspection Service, U.S. Department of Agriculture, Washington, D.C.
7. ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
8. Sadovski (1977) J. Food Technol. 12:85.
9. Edel and Kampelmacher (1973) Bull. W.H.O. 48:167.

PRESENTAZIONE	Categoria	Confezione	Ref.
Buffered Peptone Water	Terreni in provette	20 x 9 ml	24199
Buffered Peptone Water	Terreni in provette	100 x 9 ml	26199
Buffered Peptone Water	Terreni in provette	20 x 10 ml	24099
Buffered Peptone Water (Double Concentration)	Terreni in provette	20 x 9 ml	24463
Buffered Peptone Water	Terreni in flaconi	6 x 90 ml	414030
Buffered Peptone Water	Terreni in flaconi	25 x 90 ml	454030
Buffered Peptone Water	Terreni in flaconi	6 x 200 ml	412090
Buffered Peptone Water	Terreni in flaconi	6 x 225 ml	414020
Buffered Peptone Water	Terreni in flaconi	25 x 225 ml	451402
Buffered Peptone Water - Bags	Terreni in sacche	3 x 3 litri	499030
Buffered Peptone Water - Bags	Terreni in sacche	3 x 5 litri	499035
Buffered Peptone Water	Terreni disidratati	100 g	621014
Buffered Peptone Water	Terreni disidratati	500 g	611014
Buffered Peptone Water	Terreni disidratati	5 kg	6110145

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