# IODINE MKTT SOLUTION

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## DESCRIPTION

IODINE MKTT SOLUTION is a supplement for the detection of *Salmonella* spp used for enrichment of MULLER KAUFFMANN TETRATHIONATE BROTH BASE ref. 610239 or 620239.

## KIT CONTENT

Each kit contains:

- 10 tubes of IODINE MKTT SOLUTION with 10 mL of solution
- 1 instructions sheet

## COMPOSITION

### IODINE MKTT SOLUTION

	Content / tube	Content / l of medium
Potassium iodide	2.5 g	5.0 g
Iodine	2.0 g	4.0 g

## PROCEDURE OF USE

Aseptically add the content of one tube (10 ml) to 500 ml of MULLER KAUFFMANN TETRATHIONATE BROTH BASE ref. 610239 or 620239, boiled, cooled to 45-50°C and added with NOVOBIOCIN MKTT Supplement. Carefully mix. Distribute into sterile tubes.

## TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of MULLER KAUFFMANN TETRATHIONATE BROTH BASE ref. 610239 or 620239.

## QUALITY CONTROL

Microbial control.

Prepare the tubes using as base MULLER KAUFFMANN TETRATHIONATE BROTH BASE ref. 610239 or 620239 enriched with NOVOBIOCIN MKTT Supplement (1 bottle in 500 ml of medium) and IODINE MKTT SOLUTION. The tubes are seeded with the strains indicated in the microbiological control table.

Incubation conditions: 24 ± 3 h at 37 ± 1°C.

Microbial control:

Control strains		Growth
<i>Salmonella typhimurium</i>	ATCC 14028	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited
<i>Salmonella seftenberg</i>	ATCC 10384	Good

## STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

## WARNING AND PRECAUTIONS

The product is classifiable as hazardous under current legislation; it is recommended to consult the Safety Data Sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

## DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

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LIOFILCHEM® S.r.l.